

CURA F

VERSATILE AND REUSABLE

This reusable and user-friendly faceshield was specifically designed and developed for healthcare professionals. Compatible with all prescription glasses, CURA-F is perfectly suited to provide effective protection against viruses, droplets and splashes. Easily assembled, adjusted, and cleaned, it ensures long-lasting comfort, performance and reliability. CURA-F is designed, manufactured and assembled in the UK using locally sourced sustainable materials. Because protecting the planet means protecting lives in the long term.

- + **LIGHTWEIGHT DISINFECTABLE FRAME (POLYAMIDE)**
- + **EASY PEEL-OFF OUTER PROTECTIVE FILM**
- + **WASHABLE AND REUSABLE**
- + **PLATINUM LITE, ANTI-FOG COATING**



TECHNOLOGIES

PLATINUM LITE

VERSIONS



PFSCURFP02



PFSCURFP03



PFSCURFP04

CURA F

VERSION	REFERENCE	LENS MARKING	FRAME MARKING	COATING	STANDARD
 CLEAR	PFSCURFP02	CE 00000 CROWN Z87	CE 00000 Z87 D3	PLATINUM LITE	ANSI Z87.1-2015
 CLEAR	PFSCURFP03	CE 00000 CROWN Z87	CE 00000 Z87 D3	PLATINUM LITE	ANSI Z87.1-2015
 CLEAR	PFSCURFP04	CE 00000 CROWN Z87	CE 00000 Z87 D3	PLATINUM LITE	ANSI Z87.1-2015

 **WARNING:** Can expose you to Bisphenol-A (BPA), a carcinogen and reproductive toxicant. See www.P65Warnings.ca.gov.



PFSCURFP02
HEADBAND/TEMPLES/BUCKLES



PFSCURFP03
HEADBAND/TEMPLES




PFSCURFP04
HEADBAND/BUCKLES

GENERAL SPECIFICATIONS

LENS	PET	STRAP	ISOPRENE - WHITE
FRAME	PA6 - HEADBAND/TEMPLES/BUCKLES	USE	MULTI USE
WEIGHT	68 g	SUPPLIED WITH	

KNOW MORE
ABOUT
MARKINGS:



OR SCAN:

PACKAGING & EAN13 CODE

SINGLE: Size: 1.50 x 1.06 x 0.8 inch	INNER BOX:	OUTER PACK: Size: 22.83 x 15 x 7.9 inch
3660740013202	-	3660740213206 - 50 pcs
3660740013103	-	3660740213107 - 50 pcs
3660740013202	-	3660740213206 - 50 pcs
CUSTOM CODE: 9004900000		

CLEANING & STORAGE

CLEANING: The lenses should be cleaned regularly. Bollé Safety recommends only using the B-clean lens cleaning system.

MAINTENANCE: Pitted, scratched or damaged lenses must be replaced by replacing the complete product. If the frame is damaged also replace the complete product. After each use it is recommended to keep the glasses or goggles in a microfiber pouch or rigid case.

STORAGE: All products must be stored in a dry and clean environment.



MADE IN: UK

DATASHEET: 11/02/25

bolle-safety.com

Europe Headquarters
34 rue de la Soie
69100 Villeurbanne
France
Tél. : 00 33 (0) 4 78 85 23 64

U.K. & Ireland
Unit C83 - Barwell Business Park
Leatherhead Road - Chessington
Surrey - KT9 2NY - UK
Tél. : 00 44 (0)208 391 4700

USA Office
2810 Caribou Court, Suite 160
Carlsbad, CA 92010
United States of America
Tél. : 800 222 6553

Asia-Pacific Office
32-40 Fairchild Street
Heatherton VIC 3202
Australia
Tél. : 1300 138 098



**THE PPE EYEWEAR
SPECIALISTS**